




19 YEAR OLD FEMALE  
Gestational age: 12 weeks, 4 days



FAMILY HISTORY  
No known risk factors  
MaterniT 21 PLUS + SCA + ESS  
ordered, result negative  
SCA - Sex Chromosome Aneuploidies  
ESS - Enhanced Sequencing Series  
(Clinically relevant microdeletions and trisomies)




GENOME-FLEX ORDERED AT  
27 WEEKS GA  
Due to ultrasound findings of  
clubfoot and heart defect  
Positive: 32.25 Mb 11q interstitial deletion



**Normal 50 Kb trace  
(for comparison)**

Each number below the trace  
represents a chromosome, from  
1 to 22, X/Y. Note that the orange line  
stays relatively flat in a normal trace.



**Positive GENOME-Flex trace**

Note the significant downward  
deviation on the orange line for  
chromosome 11q signifying a loss of  
chromosomal material.


## Key points

- A traditional NIPT will only report on specific regions of the genome (typically T13/18/21, SCAs, X/Y), therefore missing chromosomal abnormalities located in regions not screened
- **MaterniT 21 PLUS now with GENOME-Flex** allows previously run MaterniT 21 PLUS samples to be rapidly re-sequenced using the deeper sequencing power of MaterniT GENOME to screen all 23 pairs of chromosomes
- As illustrated by this case study, using traditional NIPT screening for average risk patients may miss clinically relevant abnormalities on other chromosomes. However, a flexible NIPT like MaterniT 21 PLUS offers a new NIPT pathway should the patient become high-risk later in pregnancy due to suspected chromosomal abnormalities



MaterniT 21 PLUS is the only NIPT to offer GENOME-Flex, a new NIPT high risk pathway. MaterniT GENOME detects up to 30% more chromosomal information than other NIPTs<sup>2,3</sup>; detects chromosomal aneuploidies and copy number variants missed by traditional NIPT; thereby providing earlier awareness and more proactive pregnancy management options.

Chromosome 11 - 32.25 Mb 11q interstitial deletion



Ideogram from the MaterniT GENOME lab report with a closeup view of the impacted chromosomal trace provides a detailed view of the region of interest. The amplitude of the purple trace shows the deviation: in this case a loss on chromosome 11q.



## Case study 6 summary

- 19 year old average risk patient, no known risk factors
- MaterniT 21 PLUS (with SCA + ESS) ordered at 12 weeks GA - Negative result
- Subsequent ultrasound identified suspected clubfoot and heart defect
- **GENOME-Flex** ordered at 27 weeks GA (no additional blood draw required) - **Positive result - 32.25 Mb 11q interstitial deletion**
- A confirmatory diagnostic procedure may be recommended according to clinical practice and society guidelines, but this option may not be desired by the patient. This patient declined amniocentesis with microarray
- Baby delivered with multiple anomalies
- **GENOME-Flex was able to re-sequence this patient's original MaterniT 21 PLUS sample and identify clinically relevant abnormalities not detected by traditional NIPT (typically without a redraw and with a rapid 72 hour turnaround time)**

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

© 2019 Laboratory Corporation of America® Holdings. All rights reserved. | rep-1375-v1-0819 | L21187-0819-1