

### CLINICAL UPDATE Case study 5: three-way translocation



26 YEAR OLD FEMALE




17 WEEKS GA

Ultrasound findings: normal

Family history of translocation;
 mother, child, and sister




MATERNIT® GENOME ORDERED AT 17 WEEKS Positive: 14.8 Mb gain on 5q and a 16.5 Mb loss on 18q



# Normal 50 Kb trace (for comparison)

Each number represents a chromosome, from 1 to 22, X, Y. Note that the orange line stays relatively flat in a normal trace.



#### Positive MaterniT GENOME trace


Note the significant upward deviation on the orange line for chromosome 5q and downward deviation for chromosome 18q, signifying a gain of material on chromosome 5 and loss of material on chromosome 18.



#### Key points

- MaterniT GENOME correctly identified complex chromosomal abnormalities consistent with an unbalanced translocation, confirmed by ultrasound and diagnostic testing with microarray
- Translocation carriers have a significant risk to experience early pregnancy miscarriage as well as an increased risk of fetal abnormalities due to inheritance of an unbalanced version of the translocation.<sup>1</sup> MaterniT GENOME is a genome-wide NIPT screening option for couples who carry a balanced translocation
- As illustrated by this case study, using traditional NIPT and screening for only common aneuploidies (T13/18/21) with cfDNA may miss clinically relevant abnormalities on other chromosomes, potentially delivering false reassurance

MaterniT GENOME detects up to 30% more chromosomal information than other NIPTs<sup>2,3</sup>; detects chromosomal aneuploidies missed by traditional NIPT; thereby providing earlier awareness and more proactive pregnancy management options.



Ideograms from the MaterniT GENOME lab report with close-up views of the impacted chromosomal trace provide a detailed view of the regions of interest. The purple trace shows the deviations: a gain of chromosome 5q and a loss of chromosome 18q. (Note the purple trace in relation to the blue trace.)

## ) Case study 5 summary

- 17 weeks GA Ultrasound findings: normal, family history of translocation
- Mother carries two independent translocations, including a three-way translocation: t(5;10;18)(q33q11.2q21.3) and t(9;13)(p22q22).
  Mother also has one child who carries the three-way translocation (balanced) and a sister who has an unbalanced form of the t(9;13) translocation
- 17 weeks GA MaterniT GENOME ordered; positive for gain on chromosome 5q and a loss on chromosome 18q
- 18 weeks GA Amniocentesis with microarray confirmed 5q duplication and 18q deletion
- 20 weeks GA Subsequent ultrasound findings: increased nuchal fold, overlapping fingers, edema of the scalp. Fetal echo revealed ASD, VSD and possible anomalous pulmonary venous return
- 25 weeks GA Ultrasound findings: ventriculomegaly

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

